

NYLON FLOCKED BRUSH WITH BREAKABLE ABS HOLDER

Intended purpose:

Sterile invasive product for use through a natural orifice for collection and transport of biological samples for subsequent microbiological analysis by specialized healthcare personnel under appropriate hygienic conditions. Swabs are medical devices (MD) and are covered by CE0318. Single-use product. The sterilization method is indicated on the primary packaging.

Presentation:

The swab consists of a breakable acrylonitrile butadiene styrene (ABS) holder and a flocked nylon head shaped for nasopharyngeal sampling. Swabs may be individually packaged in a tube or peel-pack or flow-pack container. The swab is intended for use in hospital settings.

Precautions:

The use of the brush does not present any difficulty for the intended users.

Not suitable for any application other than its intended use. Do not use if there are signs of breakage or deterioration of the packaging. If the unit packaging has lost its integrity, do not use the swab as sterility cannot be guaranteed. Do not use if the swab is dirty. Discard the swab if it is broken or missing the head. Reuse of the product may cause infection. Handle with care. Although unlikely, applying too much pressure could break the swab. In case of breakage, discard immediately. Swabs do not require special storage conditions. Once used, swabs must be decontaminated or treated as contaminated hospital waste. It is the responsibility of each laboratory to handle, treat, and dispose of waste in accordance with current legislation. Unused swabs can be considered non-hazardous products and can be disposed of according to this criterion. Do not use if the expiration date has passed. Do not use if any component is missing or if the printed information on the product is not displayed correctly.

REF.	DESCRIPTION
310217	FLOCKED SWAB NASOPH. BP 80MM FLOW

Symbol glossary:

 Catalogue number	 Batch number	 Use-by-date	 Sterile by irradiation	 Unique sterile barrier system
 Unique device identifier	 Medical Device	 Do not use if package is damaged	 Do not re-use	 Manufacturing date
 Keep away from sunlight	 Manufacturer	 CE marking	 Consult instructions for use on the website www.deltalabgroup.com/eifus	 Quantity

In case of a serious incident* related to the product, notify to Deltalab, S.L. as well as the competent authority of the State in which the user is established. *A "serious incident" is understood as one that entails the death, or serious deterioration of the health of the patient or user or a serious threat to public health.

