

60ml vacuum container with sterility seal

An innovative, practical and sustainable product



Environmental friendly

 **deltalab**
member of **SCGP**



Goodbye to the single wrapping bag!

New urine sample collection container with vacuum system manufactured by Deltalab under sterile conditions (Sterile A).

Its innovative design, with a **patented molded seal**, guarantees sterility and inviolability, thus eliminating the need for a single bag.

The 60 ml volume allows to optimize the collecting space adapting to health protocols and at the same time achieving clear advantages at a functional, logistical and environmental level.

With the vacuum system, a safe transfer of urine from the container to the tube is achieved, preventing the risk of loss, splashes, or cross-contamination during the transfer.



code	description	capacity	max. vol.	case qty.
409531	60ml container with seal	60ml	82ml	650
409520	Container with seal + vacuum tube set	60ml 9ml	82ml 9ml	250



Directive 98/79/EC on *in vitro* diagnostic medical devices

The main product improved features are:

- ✓ The plastic wrapping is eliminated thanks to the **incorporation of the sterility seal in the container itself** and this implies a significant reduction in plastic waste of approximately 32%.
- ✓ Compared to 120ml containers, a **large amount of raw material is economized**.
- ✓ Thanks to its smaller size, considerable **space savings** are achieved in transport and storage.
- ✓ This new container has a **wide mouth diameter of 56mm** which offers great comfort to the user.
- ✓ The **label** is also a key element of this new product, containing all the necessary and regulatory information for its use.
- ✓ Its innovative design with an **overmolded** (not threaded) **needle** prevents any type of leak.
- ✓ The **low position of the needle** minimizes the dead volume remaining in the flask.
- ✓ It is **compatible** with the vast majority of vacuum urine tubes on the market.